



## Ca<sup>++</sup>Mg<sup>++</sup> —ATPase Activity Assay kit

---

Ca<sup>++</sup> Mg<sup>++</sup>-ATPase is widely distributed in plants, animals, microorganisms and cells, and can catalyze the hydrolysis of ATP to produce ADP and inorganic phosphorus.

According to the Ca<sup>++</sup> Mg<sup>++</sup>-ATP enzyme decomposes ATP to generate ADP and inorganic phosphorus, the level of ATPase activity is determined by measuring the amount of inorganic phosphorus.

Catalog No.	250096
Size	50 Assays / 100 Assays
Product Category	Enzymatic / Metabolic Assay
Detection Method	Visible spectrophotometry
Storage/Stability	-20°C/1 year
Shipping	Gel Packs