



Penicillin-Streptomycin-Amphotericin B [100x] (sterile)

Composition: Penicillin (10kU/mL), Streptomycin (10mg/mL), Amphotericin B (25ug/mL) in 0.9% NaCl or PBS buffer.

The 100x concentrated antibiotic solution can be added directly in a ratio of 5mL per 500mL, or 10mL per 1L of the cell culture medium; after the preparation is completed then filter and sterilized.

Catalog No.	523026
Size	100mL
Product Category	Cell Culture Reagent
Storage/Stability	-20°C/1 year
Shipping	Gel Packs

www.realgenelabs.com

For Research Use only