



Penicillin-Streptomycin-L-Glutamine [100x] (sterile)

Composition: Penicillin (10kU/mL), Streptomycin (10mg/mL) and L-glutamine (29.2 mg/mL) in 10mM citrate buffer.

The 100x concentrated antibiotic solution can be added directly in a ratio of 5mL per 500mL, or 10mL per 1L of the cell culture medium; after the preparation is completed then filter and sterilized.

Catalog No.	523036
Size	100mL
Product Category	Cell Culture Reagent
Storage/Stability	-20°C/1 year
Shipping	Gel Packs