

Bovine F2 (Coagulation Factor II) ELISA kit

This kit applies sandwich ELISA method for the quantitative detection of F2 (Coagulation Factor II) in bovine samples. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Cattle F2. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Cattle F2. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Cattle F2, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm \pm 10nm. The concentration of Cattle F2 in the samples is then determined by comparing the OD of the samples to the standard curve.

Catalog No.	3128258
Size	96-Wells
Product Category	ELISA (Quantitative)
Reactivity	Bovine
Sample	Serum
Assay Method	Sandwich ELISA
Assay Duration	3.5 hours
Sensitivity	0.288 ng/mL
Standard Curve Range	0.63-40 ng/mL
Storage/Stability	-20°C/1 year; 4°C/6 months

www.realgenelabs.com

For Research Use only