

Bovine HSPA1A (Heat Shock 70kDa Protein 1A) ELISA kit

This kit applies sandwich ELISA method for the quantitative detection of HSPA1A (Heat Shock 70kDa Protein 1A) in bovine samples. The microtiter plate provided in this kit has been precoated with an antibody specific to Cattle HSPA1A. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Cattle HSPA1A. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Cattle HSPA1A, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ±10nm. The concentration of Cattle HSPA1A in the samples is then determined by comparing the OD of the samples to the standard curve.

| Catalog No. | 3128002 |
|----------------------|--|
| Size | 96-Wells |
| Product Category | ELISA (Quantitative) |
| Reactivity | Bovine |
| Sample | Serum; Plasma; Tissue Homogenates; Cell Lysates; Cell Culture Supernates; Other Biological Fluids |
| Assay Method | Sandwich ELISA |
| Assay Duration | 3.5 hours |
| Sensitivity | 0.59 ng/mL |
| Standard Curve Range | 1.57-100 ng/mL |
| Storage/Stability | -20°C/1 year; 4°C/6 months |