

## Bovine ITaC (Interferon Inducible T-Cell alpha Chemoattractant) ELISA kit

---

This kit applies sandwich ELISA method for the quantitative detection of ITaC (Interferon Inducible T-Cell alpha Chemoattractant) in bovine samples. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Cattle ITaC. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Cattle ITaC. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Cattle ITaC, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm  $\pm$ 10nm. The concentration of Cattle ITaC in the samples is then determined by comparing the OD of the samples to the standard curve.

Catalog No.	3127977
Size	96-Wells
Product Category	ELISA (Quantitative)
Reactivity	Bovine
Sample	Serum; Plasma; Other Biological Fluids
Assay Method	Sandwich ELISA
Assay Duration	3.5 hours
Sensitivity	15.46 pg/mL
Standard Curve Range	46.88-3000 pg/mL
Storage/Stability	-20°C/1 year; 4°C/6 months