

## Bovine PPAR-gamma (Peroxisome Proliferator Activated Receptor gamma) ELISA kit

---

This kit applies sandwich ELISA method for the quantitative detection of PPAR-gamma (Peroxisome Proliferator Activated Receptor gamma) in bovine samples. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Cattle PPAR-gamma. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Cattle PPAR-gamma. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Cattle PPAR-gamma, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm  $\pm$ 10nm. The concentration of Cattle PPAR-gamma in the samples is then determined by comparing the OD of the samples to the standard curve.

Catalog No.	3126488
Size	96-Wells
Product Category	ELISA (Quantitative)
Reactivity	Bovine
Sample	Serum; Plasma; Tissue Homogenates; Cell Lysates; Cell Culture Supernates; Other Biological Fluids
Assay Method	Sandwich ELISA
Assay Duration	3.5 hours
Sensitivity	0.133 ng/mL
Standard Curve Range	0.32-20 ng/mL
Storage/Stability	-20°C/1 year; 4°C/6 months

