

Bovine TSP-1 (Thrombospondin-1) ELISA kit

This kit applies sandwich ELISA method for the quantitative detection of TSP-1 (Thrombospondin-1) in bovine samples. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Cattle TSP-1. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Cattle TSP-1. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Cattle TSP-1, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ±10nm. The concentration of Cattle TSP-1 in the samples is then determined by comparing the OD of the samples to the standard curve.

| Catalog No. | 3125945 |
|----------------------|----------------------------|
| Size | 96-Wells |
| Product Category | ELISA (Quantitative) |
| Reactivity | Bovine |
| Sample | Serum; Plasma |
| Assay Method | Sandwich ELISA |
| Assay Duration | 3.5 hours |
| Sensitivity | 30 pg/mL |
| Standard Curve Range | 78.13-5000 pg/mL |
| Storage/Stability | -20°C/1 year; 4°C/6 months |