

## Bovine uPA (Plasminogen Activator, Urokinase) ELISA kit

---

This kit applies sandwich ELISA method for the quantitative detection of uPA (Plasminogen Activator, Urokinase) in bovine samples. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Cattle uPA. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Cattle uPA. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Cattle uPA, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm  $\pm$ 10nm. The concentration of Cattle uPA in the samples is then determined by comparing the OD of the samples to the standard curve.

Catalog No.	3125864
Size	96-Wells
Product Category	ELISA (Quantitative)
Reactivity	Bovine
Sample	Serum; Plasma; Other Biological Fluids
Assay Method	Sandwich ELISA
Assay Duration	3.5 hours
Sensitivity	0.061 ng/mL
Standard Curve Range	0.16-10 ng/mL
Storage/Stability	-20°C/1 year; 4°C/6 months