

Canine HMGB-1 (HighMobility Group Protein B1) ELISA kit

This kit applies sandwich ELISA method for the quantitative detection of HMGB-1 (HighMobility Group Protein B1) in canine samples. The microtiter plate provided in this kit has been precoated with an antibody specific to Dog HMGB-1. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Dog HMGB-1. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Dog HMGB-1, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm ±10nm. The concentration of Dog HMGB-1 in the samples is then determined by comparing the OD of the samples to the standard curve.

| Catalog No. | 3145992 |
|----------------------|--|
| Size | 96-Wells |
| Product Category | ELISA (Quantitative) |
| Reactivity | Canine |
| Sample | Serum; Plasma; Other Biological Fluids |
| Assay Method | Sandwich ELISA |
| Assay Duration | 3.5 hours |
| Sensitivity | 2.73 pg/mL |
| Standard Curve Range | 6.25-400 pg/mL |
| Storage/Stability | -20°C/1 year; 4°C/6 months |