

Canine SPA2 (Surfactant Associated Protein A2) ELISA kit

This kit applies sandwich ELISA method for the quantitative detection of SPA2 (Surfactant Associated Protein A2) in canine samples. The microtiter plate provided in this kit has been precoated with an antibody specific to Dog SPA2. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Dog SPA2. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Dog SPA2, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm \pm 10nm. The concentration of Dog SPA2 in the samples is then determined by comparing the OD of the samples to the standard curve.

Catalog No.	3149772
Size	96-Wells
Product Category	ELISA (Quantitative)
Reactivity	Canine
Sample	Tissue Homogenates; Cell Lysates; Cell Culture Supernates; Other Biological Fluids
Assay Method	Sandwich ELISA
Assay Duration	3.5 hours
Sensitivity	0.062 ng/mL
Standard Curve Range	0.16-10 ng/mL
Storage/Stability	-20°C/1 year; 4°C/6 months