

Rat Pro-ANP (Pro Atrial Natriuretic Peptide) ELISA kit

This kit applies sandwich ELISA method for the quantitative detection of Pro-ANP (Pro Atrial Natriuretic Peptide) in rat samples. The microtiter plate provided in this kit has been pre-coated with an antibody specific to Rat Pro-ANP. Standards or samples are added to the appropriate microtiter plate wells then with a biotin-conjugated antibody specific to Rat Pro-ANP. Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Rat Pro-ANP, biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm \pm 10nm. The concentration of Rat Pro-ANP in the samples is then determined by comparing the OD of the samples to the standard curve.

| | |
|----------------------|---|
| Catalog No. | 3038360 |
| Size | 96-Wells |
| Product Category | ELISA (Quantitative) |
| Reactivity | Rat |
| Sample | Serum; Plasma; Tissue Homogenates; Cell Lysates; Cell Culture Supernates; Other Biological Fluids |
| Assay Method | Sandwich ELISA |
| Assay Duration | 3.5 hours |
| Sensitivity | 15.6 pg/mL |
| Standard Curve Range | 62.5-4000 pg/mL |
| Storage/Stability | -20°C/1 year; 4°C/6 months |

