



50bp DNA Ladder [250uL]

A ready-to-use 50bp DNA Ladder consists of 14 DNA fragments ranging from 50bp-1000bp. Recommended to load 5ul each lane; optimized to create a sharp and clear band with the enhanced intensity.

Highly stable at room temperature. If kept in the freezer, multiple freeze-thaw cycles will not affect the quality of DNA markers.

Catalog No.	530020
Size	100 Loads
Product Category	Nucleic Acid Electrophoresis
Storage/Stability	-20°C/2 years
Shipping	Gel Packs