



## 100bp DNA Ladder [250uL]

---

A ready-to-use 100bp DNA Ladder consists of 10 DNA fragments ranging from 100bp-1000bp. Recommended to load 5ul each lane; optimized to create a sharp and clear band with the enhanced intensity.

Highly stable at room temperature. If kept in the freezer, multiple freeze-thaw cycles will not affect the quality of DNA markers.

|                   |                              |
|-------------------|------------------------------|
| Catalog No.       | 530121                       |
| Size              | 100 Loads                    |
| Product Category  | Nucleic Acid Electrophoresis |
| Storage/Stability | -20°C/2 years                |
| Shipping          | Gel Packs                    |