



## 100bp PLUS DNA Ladder [250uL]

---

A ready-to-use 100bp PLUS DNA Ladder consists of 14 DNA fragments ranging from 100bp-3000bp. Recommended to load 5ul each lane; optimized to create a sharp and clear band with the enhanced intensity.

Highly stable at room temperature. If kept in the freezer, multiple freeze-thaw cycles will not affect the quality of DNA markers.

Catalog No.	530021
Size	100 Loads
Product Category	Nucleic Acid Electrophoresis
Storage/Stability	-20°C/2 years
Shipping	Gel Packs

[www.realgenelabs.com](http://www.realgenelabs.com)

For Research Use only