

Acetyl CoA (AcCoA) Rapid ELISA kit

This Rapid ELISA kit applies Sandwich enzyme immunoassay technique for the quantitative detection. The micro ELISA plate provided in this kit has been pre-coated with an antibody specific to Acetyl CoA (AcCoA) . Standards or samples are added to the appropriate micro ELISA plate wells then with a biotin-conjugated antibody specific to Acetyl CoA (AcCoA) . Next, Avidin conjugated to Horseradish Peroxidase (HRP) is added to each microplate well and incubated. After TMB substrate solution is added, only those wells that contain Acetyl CoA (AcCoA) , biotin-conjugated antibody and enzyme-conjugated Avidin will exhibit a change in color. The enzyme-substrate reaction is terminated by the addition of sulphuric acid solution and the color change is measured spectrophotometrically at a wavelength of 450nm \pm 10nm. The concentration of Acetyl CoA (AcCoA) in the samples is then determined by comparing the OD of the samples to the standard curve.

Catalog No.	3510014
Size	96-Well Tests
Product Category	Rapid ELISA
Reactivity	Various
Sample	Celll Lysate; Tissue Lysate
Assay Method	Sandwich ELISA
Assay Duration	1.5H
Sensitivity	0.119 nmol/mL
Standard Curve Range	0.32-20 nmol/mL
Storage/Stability	-20°C/1 year; 4°C/6 months