

## Mouse Vascular Endothelial Growth Factor Receptor 2 (VEGFR-2/KDR) Rapid ELISA kit

This Rapid ELISA kit applies Sandwich enzyme immunoassay technique for the quantitative detection. The micro ELISA plate provided in this kit has been pre-coated with an antibody specific to Mouse VEGFR-2/KDR. Standards or samples are added to the micro ELISA plate wells and combined with the specific antibody. Then a biotinylated detection antibody specific for Mouse VEGFR-2/KDR and Avidin-Horseradish Peroxidase (HRP) conjugate are added successively to each micro plate well and incubated. Free components are washed away. The substrate solution is added to each well. Only those wells that contain Mouse VEGFR-2/KDR, biotinylated detection antibody and Avidin-HRP conjugate will appear blue in color. The enzyme-substrate reaction is terminated by the addition of stop solution and the color turns yellow. The optical density (OD) is measured spectrophotometrically at a wavelength of 450 nm ± 2 nm. The OD value is proportional to the concentration of Mouse VEGFR-2/KDR. You can calculate the concentration of Mouse VEGFR-2/KDR in the samples by comparing the OD of the samples to the standard curve.

Catalog No.	3531933
Size	96-Well Tests
Product Category	Rapid ELISA
Reactivity	Mouse
Sample	Serum; Plasma; Tissue Homogenate; Other Biological Samples
Assay Method	Sandwich ELISA
Assay Duration	1.5H
Sensitivity	0.1 ng/mL
Standard Curve Range	0.16-10 ng/mL

C+	/C+-	La Hilland
Stora	ge/Sta	IDIIITY

-20°C/1 year; 4°C/6 months

www.realgenelabs.com

For Research Use only