

Rat Peroxisome Proliferator Activated Receptor gamma (PPAR-gamma) Rapid ELISA kit

This Rapid ELISA kit applies Sandwich enzyme immunoassay technique for the quantitative detection. The micro ELISA plate provided in this kit has been pre-coated with an antibody specific to Rat PPAR-gamma. Standards or samples are added to the micro ELISA plate wells and combined with the specific antibody. Then a biotinylated detection antibody specific for Rat PPAR-gamma and Avidin-Horseradish Peroxidase (HRP) conjugate are added successively to each micro plate well and incubated. Free components are washed away. The substrate solution is added to each well. Only those wells that contain Rat PPAR-gamma, biotinylated detection antibody and Avidin-HRP conjugate will appear blue in color. The enzyme-substrate reaction is terminated by the addition of stop solution and the color turns yellow. The optical density (OD) is measured spectrophotometrically at a wavelength of 450 nm \pm 2 nm. The OD value is proportional to the concentration of Rat PPAR-gamma. You can calculate the concentration of Rat PPAR-gamma in the samples by comparing the OD of the samples to the standard curve.

| Catalog No. | 3541782 |
|----------------------|---|
| Size | 96-Well Tests |
| Product Category | Rapid ELISA |
| Reactivity | Rat |
| Sample | Tissue Homogenate; Celll Lysate; Cell culture supernatant |
| Assay Method | Sandwich ELISA |
| Assay Duration | 1.5H |
| Sensitivity | 0.1 ng/mL |
| Standard Curve Range | 0.16-10 ng/mL |
| Storage/Stability | -20°C/1 year; 4°C/6 months |