



REAL-Detect Hazelnut Oleosin gene qPCR kit

REAL-Detect Hazelnut Oleosin gene qPCR kit applies qualitative fluorescence PCR method for gene detection. The kit is supplied with reaction mix, positive control, negative control and re-dissolved diluent. The components supplied in the kit are in lyophilized format (stable at 25-30°C for 2 weeks).

Kit is compatible with all the standard Real-Time PCR systems, including:

ABI: 7000, 7300, 7500, 7900, Step One, Step One Plus, QuantStudio

Agilent: Mx3000P, 3005P

Bio-Rad: CFX 96, iCycler iQ4, iCycleriQ5

Cepheid: Smart Cyclo-II

MJ: Option2, Chromo4

Qiagen: Rotor Gene 6000

Roche: LightCycler680

| | |
|-------------------|----------------------|
| Catalog No. | 950003 |
| Size | 48 Tests |
| Product Category | Transgenic Plant PCR |
| Storage/Stability | -20°C / 1 year |
| Shipping | Ambient |