



REAL-Plex Human Papillomavirus 16 & 18 (HPV 16/18) qPCR Panel

REAL-Plex Human Papillomavirus 16 & 18 (HPV 16/18) qPCR Panel applies qualitative multiplex fluorescence PCR method for pathogen detection. The kit is supplied with reaction mix, positive control, negative control and re-dissolved diluent. The components supplied in the kit are in lyophilized format (stable at 25-30°C for 2 weeks).

Kit is compatible with all the standard Real-Time PCR systems, including:

ABI: 7000, 7300, 7500, 7900, Step One, Step One Plus, QuantStudio

Agilent: Mx3000P, 3005P

Bio-Rad: CFX 96, iCycler iQ4, iCycleriQ5

Cepheid: Smart Cyclor-II

MJ: Option2, Chromo4

Qiagen: Rotor Gene 6000

Roche: LightCycler680

Catalog No.	991012
Size	
Product Category	Multiplex PCR
Storage/Stability	-20°C / 1 year
Shipping	Gel Packs
Pack Size / Price	Inquire