

Nissl stain kit (tar violet method)

NissI body is a triangular or oval small lumps distributed in the cytoplasm of nerve cells, which can be stained purple-blue by basic dyes. All kinds of nerve cells contain NissI bodies, but their shape, quantity and distribution position are often different. NissI bodies are also present in dendrites, but not in axons and axon mounds of inclusions. NissI body is an important part of protein synthesis in neurons. When neurons are stimulated, the NissI body in inclusions will be significantly reduced. Changes in NissI bodies can be well visualized by tar violet staining. This product is easy to operate, stable in staining, and has a wide range of applications. It can be used for staining NissI bodies and neurons in paraffin tissue sections.

Catalog No.	260028
Size	2 x 100mL
Product Category	Histochemical Stain
Storage/Stability	Ambient/6 months
Shipping	Ambient

www.realgenelabs.com

For Research Use only