



## Native Lysis Buffer

---

Native lysis buffer contains a non-ionic detergent that lyses the cells and releases cytosolic proteins and soluble membrane proteins, nuclear proteins under non-denaturing conditions. The lytic protein product under non-denaturing conditions maximizes the properties and functions of the protein, such as antigen-antibody binding or enzymatic activity, and is therefore suitable for co-immunoprecipitation. In addition, some antibodies have a higher binding capacity for non-denatured proteins or can only recognize antigenic sites of non-denatured proteins, in which case non-denaturing cleavage should be used.

### Kit Components:

1. Native Lysis Buffer: 100mL
2. PMSF (100mM): 1.5mL

Catalog No.	501001
Size	100mL
Product Category	Protein Extraction
Storage/Stability	2 ~ 8°C/1 year
Shipping	Gel Packs