



CRISPR-Cas13a DNA Lateral Flow Detection Kit (2-step)

Cas13a forms a functional complex with crRNA, "gliding" over the target nucleic acid sequence. Upon successful pairing, Cas13a is specifically activated and trans-cleaves the surrounding reporter molecule. This kit combines isothermal amplification technology with Cas13a, boasting high sensitivity, specificity, and a high signal-to-noise ratio. It has been widely used in molecular diagnostics, enabling rapid and accurate detection of pathogens.

This kit applied isothermal amplification & CRISPR/Cas13a, enables DNA detection in two steps, and the results can be observed visually using a chromatographic test strip.

Catalog No.	519226
Size	
Product Category	Isothermal Amplification
Pack Size	96 rxns
Storage/Stability	-20°C/1 year
Shipping	Gel Packs