



Lateral Flow Strip (for RPA)

This product detects FITC/FAM and biotin dual-labeled gene amplification products based on lateral flow technology using gold particles. One primer/probe labeled with biotin and the other with fluorescein isothiocyanate (FITC) or 6-carboxyfluorescein (6-FAM), ensuring that both labels can be incorporated into the double-stranded amplification product at the same time. Recommend to use with EZassay™ products Isothermal DNA amplification Nfo kit.

Catalog No.	519253
Size	
Product Category	Isothermal Amplification
Pack Size	48 rxns / 96 rxns
Storage/Stability	2 ~ 8°C/1 year (Protect from moisture and direct sunlight)
Shipping	Ambient