



100bp DNA Ladder

A ready-to-use 100bp DNA is provided with 1x Loading Buffer for tracking of sample DNA migration, contains 10 fragments: 100bp, 200bp, 300bp, 400bp, 500bp, 600bp, 700bp, 800bp, 900bp and 1000bp. 5uL of this product contains ~50ng of each band, except the 500bp fragment which contains ~120ng.

Recommended to load 5ul each lane; optimized to create a sharp and clear band with the enhanced intensity.

Highly stable at room temperature. If kept in the freezer, multiple freeze-thaw cycles will not affect the quality of DNA marker.

Catalog No.	530121
Size	250uL / 500uL
Product Category	DNA Ladder
Storage/Stability	-20°C/2 years
Shipping	Gel Packs