



## Coomassie Blue Staining/Destaining Kit

---

Coomassie Blue Staining/Destaining Kit uses coomassie Brilliant Blue staining and destaining method. It uses Coomassie Brilliant Blue G250 as the dye and can be used for protein electrophoresis gels such as SDS-PAGE or non-denaturing PAGE. Routine detection or staining and destaining of residual protein on PAGE gel after Western transfer.

### Kit Components:

1. Coomassie Brilliant Blue Staining Solution: 100mL
2. Coomassie Brilliant Blue Destaining Solution: 500mL

Catalog No.	502202
Size	Kit / WSP
Product Category	SDS-PAGE Buffer & Stain
Storage/Stability	Ambient/1 year
Shipping	Ambient