



## Silver Stain kit

---

The silver staining of protein bands is based on the combination of various groups in the protein (such as mercapto groups, carbon groups, etc.) Ribbon stained dark brown. PAGE Gel Silver Stain kit is mainly used for polyacrylamide gel electrophoresis staining. Its sensitivity is higher than that of Coomassie Brilliant Blue. The whole process of this product is simple to operate and has high sensitivity. It can detect protein bands below 10ng.

### Kit Components:

1. Staining solution A: 200mL;
2. Staining solution B: 1mL;
3. Staining solution C: 100mL;
4. Chromogenic solution D: 100mL;
5. Reagent E: 2g

Catalog No.	502204
Size	Kit / WSP
Product Category	SDS-PAGE Buffer & Stain
Storage/Stability	Ambient/1 year
Shipping	Ambient