



1x Tris-Tricine-SDS Buffer

Tris-Tricine Gel Buffer separates proteins or peptides with molecular weights of 110kD by using SDS-PAGE. The Tris-Tricine is ready-to-use gel buffer consist all the necessary reagents for the preparation of TricineSDS-PAGE gel.

Catalog No.	503016
Size	500mL / Size1 / WSP
Product Category	SDS-PAGE Buffer & Stain
Storage/Stability	2 ~ 8°C/1 year
Shipping	Ambient

www.realgenelabs.com

For Research Use only