



## Rapid PAGE-Gel Preparation kit [6%]

---

Rapid PAGE-Gel Preparation kit contains various reagents for PAGE Gel preparation, only the experimental instruments and ddH<sub>2</sub>O are needed to prepare. There is no need to calculate the required amount of solution, no need to dilute, no need to add TEMED in the process of gel preparation, just add modified coagulant to gel. The red dye makes the wells of the stacking gel clearer and easier to load prepared samples. This kit is enough for 30-50 pieces of normal size PAGE gel. The kit is simple to operate with simple in gel preparation, fast and safe in gel distribution.

### Kit Components:

1. SDS-PAGE Separating Gel Premix;
2. SDS-PAGE Stacking Gel Premix;
3. SDS-PAGE gel coagulant;
4. SDS-PAGE gel accelerator

|                   |                                                 |
|-------------------|-------------------------------------------------|
| Catalog No.       | 503207                                          |
| Size              | 25 preps / 50 preps / Size1 / WSP / Size2 / WSP |
| Product Category  | SDS-PAGE Buffer & Stain                         |
| Storage/Stability | 2 ~ 8°C/1 year                                  |
| Shipping          | Gel Packs                                       |