

Glutamic-Oxalacetic Transaminase/Aspartate Aminotransferase Activity Assay kit

GOT/AST (2.6.1.1) widely exists in animals, plants, microorganisms and cultured cells. It catalyzes the reversible amino reaction and is an important enzyme for amino acid metabolism. In addition, GOT has the highest content in myocardial cells, and is generally used as an auxiliary examination for myocardial infarction and myocarditis in clinical practice. The serum concentration can also increase when the liver is damaged.

GOT catalyzes the transamination reaction of alpha-homoglutarate and aspartic acid to produce glutamic acid and oxaloacetic acid. Oxaloacetic acid further decarboxylates itself to produce pyruvate; pyruvate can be combined with 2,4-dinitro 2,4-Dinitrophenylhydrazine reacts to produce 2,4-dinitrophenylhydrazone, which appears brown-red under alkaline conditions; the GOT enzyme activity can be calculated by measuring the change in absorbance at 505nm.

Catalog No.	250156
Size	50 Assays / 100 Assays
Product Category	Colorimetric Assay
Detection Method	Visible spectrophotometry
Storage/Stability	2 ~ 8°C/1 year
Shipping	Gel Packs