



## Human AGTR1 (Angiotensin II Receptor 1) ELISA Kit

|                   |                                                                                                                           |
|-------------------|---------------------------------------------------------------------------------------------------------------------------|
| Catalog No.       | 3021250                                                                                                                   |
| Size              | 96 Tests                                                                                                                  |
| Product Category  | ELISA (Quantitative)                                                                                                      |
| Reactivity        | Human                                                                                                                     |
| Sample Type       | Tissue homogenates, Cell lysates & other biological fluids                                                                |
| Assay Method      | Sandwich                                                                                                                  |
| Assay Duration    | 3.5 hours                                                                                                                 |
| Sensitivity       | 0.057 ng/mL                                                                                                               |
| Standard          | 10 ng/mL                                                                                                                  |
| Detection Range   | 0.16-10 ng/mL                                                                                                             |
| Synonyms          | AT1R; AG2S; AG2-S; AGTR1A; AGTR1B; AT1; AT1B; AT2R1; AT2-R1; AT2R1A; AT2R1B; HAT1R; AT1 Receptor; Angiotensin Receptor 1B |
| Storage/Stability | -20°C/1 year; 4°C/6 months                                                                                                |