



Protein G Resin column [5mL]

The 5mL cartridges are pre-filled with Protein G Resin FF (Fast Flow), a highly cross-linked 4% agarose gel matrix that allows the purification of monoclonal and polyclonal antibodies at relatively high flow rates.

Catalog No.	508052
Size	Ea
Product Category	Afiinity Purification
Matrix	Highly Cross-Linked 4% Agarose Beads
Ligand	Recombinant Protein G
Binding Capacity	? 20mg Human IgG / mL resin
Bead Size	90-150um
pH Stability	3-10
Storage Buffer	1x PBS containing 20% Ethanol
Storage/Stability	2-8°C/2 years
Shipping	Ambient