



## Ni Agarose Resin column [5mL]

---

The 5mL cartridges are pre-filled with Ni-NTA Resin FF (Fast Flow), a highly cross-linked 6% agarose gel matrix that allows the purification of Histidine-Tagged (6x His-Tagged) proteins

Catalog No.	508054
Size	Ea
Product Category	Afiinity Purification
Matrix	Highly Cross-Linked 6% Agarose
Ligand	Nitrilotriacetic Acid (NTA)
Binding Capacity	? 40mg 6x His-Tagged Protein / mL resin
Bead Size	90-150um
pH Stability	2-14
Storage Buffer	1x PBS containing 20% Ethanol
Storage/Stability	2-8°C/2 years
Shipping	Ambient