



## Glutathione Resin column [1mL]

---

The 1mL cartridges are pre-filled with Glutathion Resin FF (Fast Flow), a highly cross-linked 4% agarose gel matrix that allows the purification of GST-Tagged proteins.

Catalog No.	508015
Size	Ea
Product Category	Afiinity Purification
Matrix	Highly Cross-Linked 4% Agarose
Ligand	Glutathione
Binding Capacity	? 10mg GST-Tagged Protein / mL resin
Bead Size	90-150um
pH Stability	2-14
Storage Buffer	1x PBS containing 20% Ethanol
Storage/Stability	2-8°C/2 years
Shipping	Ambient