



## Penicillin-Streptomycin-Amphotericin B [100x]

---

Composition: Penicillin (10kU/mL), Streptomycin (10mg/mL), Amphotericin B (25ug/mL) in 0.9% NaCl or PBS buffer.

The 100x concentrated antibiotic solution can be added directly in a ratio of 5mL per 500mL, or 10mL per 1L of the cell culture medium; after the preparation is completed then filter and sterilized.

|                   |                      |
|-------------------|----------------------|
| Catalog No.       | 523032               |
| Size              | 100mL                |
| Product Category  | Cell Culture Reagent |
| Storage/Stability | -20°C/1 year         |
| Shipping          | Gel Packs            |

[www.realgenelabs.com](http://www.realgenelabs.com)

For Research Use only