



Ca⁺⁺Mg⁺⁺ —ATPase Activity Assay kit

Ca⁺⁺ Mg⁺⁺-ATPase is widely distributed in plants, animals, microorganisms and cells, and can catalyze the hydrolysis of ATP to produce ADP and inorganic phosphorus.

According to the Ca⁺⁺ Mg⁺⁺-ATP enzyme decomposes ATP to generate ADP and inorganic phosphorus, the level of ATPase activity is determined by measuring the amount of inorganic phosphorus.

Catalog No.	250096
Size	50 Assays / 100 Assays
Product Category	Colorimetric Assay
Detection Method	Spectrophotometry / Micro-Plate Reader
Storage/Stability	-20°C/1 year
Shipping	Gel Packs