



## 2x RNA Loading Buffer (Denatured)

---

A ready-to-use and room temperature stable buffer, to be used for loading RNA samples onto denaturing agarose gel. The buffer contains SDS, EDTA, Bromphenol Blue, Xylene Cyanol FF and Deionized formamide in recommended concentrations.

Catalog No.	530012
Size	1mL
Product Category	Nucleic Acid Electrophoresis
Storage/Stability	2 ~ 8°C/1 year
Shipping	Gel Packs

[www.realgenelabs.com](http://www.realgenelabs.com)

For Research Use only