

RNase Away (spray bottle)

RNase Away is a RNase inactivator, contains a variety of components that can effectively inactivate RNase contamination on the solid surface and ensures the cleanliness of the RNA working environment.

Efficient: completely inactivates RNase within 5 minutes on solid surfaces, with much higher potency and speed than commonly used DEPC.

Non-toxic: unlike DEPC, which is a strong carcinogen, this product does not have any toxicity on human body.

Can inactivate RNase T1, RNase H, BAL31, S1, Mung bean nuclease and other RNases and DNases.

Catalog No.	156002
Size	500mL
Product Category	Stock Solution / Buffer
Storage/Stability	Ambient/1 year
Shipping	Ambient

www.realgenelabs.com

For Research Use only