



## ELISA Stop Solution (blocking reagent for TMB substrate)

---

Ready-to-use working solution for TMB Substrate; it is used to stop TMB substrate reaction for ELISA and change the absorbance from blue (620 nm) to yellow (450 nm).

IFU: After adding TMB to the ELISA plate, add 50ul of ELISA stop solution to stop the reaction, and measure the absorbance at 450nm within 30 minutes.

Catalog No.	506018
Size	100mL / 500mL
Product Category	Immuno Reagent
Storage/Stability	Ambient/1 year
Shipping	Ambient