



Transhydrogenase-2 Activity Assay kit

TH is located on the inner membrane of mitochondria, also known as mitochondrial complex VI. It catalyzes the mutual conversion of $\text{NADH} + \text{NADP}^+$ and $\text{NAD}^+ + \text{NADPH}$, and regulates the balance of mitochondrial NAD(H) and NADP(H) . The reverse reaction is called TH-2, which catalyzes NADPH and NAD^+ to produce NADP^+ and NADH .

NADH and NADPH both have characteristic absorption at 340nm, so the hydrogen transfer reaction catalyzed by TH cannot cause a change in absorbance at 340nm. Using the synthetic substrate 3-acetylpyridine adenine dinucleotide (APAD^+) instead of NAD^+ , TH-2 catalyzes the reduction of APAD^+ to APADH . APADH has a characteristic light absorption at 375nm. The increase rate of light absorption at 375nm is measured to calculate TH- 2 Activity.

Catalog No.	250326
Size	50 Assays / 100 Assays
Product Category	Colorimetric Assay
Storage/Stability	-20°C/1 year
Shipping	Gel Packs