



## RT-LAMP Isothermal Amplification kit (w/ UDG)

---

This product aims to simplify RT-LAMP assay design with robust reaction and high tolerance of inhibitors. It can be combined with fluorescent probes or dyes, lateral flow and agarose gel detections to observe results. Non-specific products are common in isothermal amplifications including RT-LAMP. AapCas12b protein or BrCas12b protein is recommended in RT-LAMP CRISPR to increase specificity. RNaseH2 enzyme is also highly recommended to increase the specificity of RT-LAMP. It performs well in SNP detection.

This kit aims to simplify RT-LAMP assay design with robust reaction and high tolerance of inhibitors. It can be combined with fluorescent probes or dyes, lateral flow and agarose gel detections to observe results.

The kit contains uracil DNA glycosylase (UDG) to prevent carryover contamination and aerosol contamination.

Catalog No.	518311
Size	96 rxns
Product Category	Isothermal Amplification
Storage/Stability	-20°C/1 year
Shipping	Gel Packs