



RT-LAMP Isothermal Amplification kit (w/ Fluorescent Dye & UDG)

This kit contains a EvaGreen,fluorescent dye, for more convenient and rapid RNA isothermal amplifications.

The kit contains uracil DNA glycosylase (UDG) to prevent carryover contamination and aerosol contamination.

Catalog No.	518312
Size	96 rxns
Product Category	Isothermal Amplification
Storage/Stability	-20°C/1 year
Shipping	Gel Packs

www.realgenelabs.com

For Research Use only