



## RPA Isothermal DNA Amplification basic kit

---

This kit provides the reagents required for isothermal amplification of nucleic acids, which has the advantages of high sensitivity, strong specificity, and short reaction time. The reaction is carried out at a constant temperature (37 ~ 42 °C), and the amplification is usually completed within 30 minutes.

This kit completes DNA amplification only. Amplicons can be observed in agarose gel. This kit also works with CRISPR SHERLOCK and DETECTR.

|                   |                          |
|-------------------|--------------------------|
| Catalog No.       | 518401                   |
| Size              | 48 rxns / 96 rxns        |
| Product Category  | Isothermal Amplification |
| Storage/Stability | -20°C/1 year             |
| Shipping          | Gel Packs                |