



CRISPR-Cas12a RNA Fluorescent Detection kit (2-step)

Cas13a forms a functional complex with crRNA, "gliding" over the target nucleic acid sequence. Upon successful pairing, Cas13a is specifically activated and trans-cleaves the surrounding reporter molecule. This kit combines isothermal amplification technology with Cas13a, boasting high sensitivity, specificity, and a high signal-to-noise ratio. It has been widely used in molecular diagnostics, enabling rapid and accurate detection of pathogens.

This kit applied isothermal amplification & CRISPR/Cas13a, enables DNA detection in two step fluorescent signal.

Catalog No.	518712
Size	96 rxns
Product Category	Isothermal Amplification
Storage/Stability	-20°C/1 year
Shipping	Gel Packs