



CRISPR-Cas13a DNA Lateral Flow Detection kit (2-step)

Cas12a forms a functional complex with crRNA, "gliding" over the target nucleic acid sequence. Upon successful pairing, Cas12a is specifically activated and trans-cleaves the surrounding reporter molecule. This kit combines isothermal amplification technology with Cas12a, boasting high sensitivity, specificity, and a high signal-to-noise ratio. It has been widely used in molecular diagnostics, enabling rapid and accurate detection of pathogens.

This kit applied isothermal amplification & CRISPR/Cas12a, enables RNA detection in two steps, and the results can be observed visually using a chromatographic test strip.

Catalog No.	518803
Size	96 rxns
Product Category	Isothermal Amplification
Storage/Stability	-20°C/1 year
Shipping	Gel Packs

www.realgenelabs.com

For Research Use only