



## Mycoplasma Removal Reagent [1000x]

---

In the process of cell cultures, regular detection can effectively avoid mycoplasma infection, but if the cells have been infected with mycoplasma, they need to be treated in time. The best treatment is to discard the cells after high-pressure sterilization to avoid contaminating other clean cell lines. If contaminated cells are precious, mycoplasma infection must be removed. Because mycoplasma lacks the cell wall, conventional antibiotics are generally ineffective. This product is an improved mixed preparation specially used for mycoplasma elimination. It can achieve a good mycoplasma elimination effect by inhibiting the synthesis of related proteins and is not harmful to cells. Minimize the loss of your precious cells due to mycoplasma infection.

This product is supplied in 1000x contraction, which needs to be diluted to the concentration of 1x prior to use.

Catalog No.	400101
Size	1mL / 5 x 1mL
Product Category	Mycoplasma PCR