



Diff-Quik staining Kit

Diff-Quik staining is a rapid staining method used in cytology. The staining kit consists of a Fixative Agent (methanol, blue), Reagent-I (eosinophilic, orange) and Reagent-II (basophilic, blue). Generally, slides are dipped sequentially into each solution for 6 times (or left for 10-15 seconds in each solution), followed by a water rinse and drying. The required time is extremely short, and the dyeing can be completed within 90s.

Diff-Quick staining is mainly used for blood cell smear, bone marrow smear, vaginal secretion smear, exfoliated cell smear. Diff-Quick staining is ideal for batch dip dyeing with clear background and no sediment.

Kit Components:

1. Diff-Quik Fixative (100mL / 500mL)
2. Diff-Quik Reagent-I (100mL / 500mL)
3. Diff-Quik Reagent-II (100mL / 500mL)

Catalog No.	260004
Size	3 x 100mL / 3 x 500mL
Product Category	Histochemical Stain
Storage/Stability	Ambient/1 year
Shipping	Ambient