



Enhance DAB substrate kit [20x]

The kit contains Diaminobenzidine (DAB) substrate and Peroxidase (HRP) enzyme for the immunostaining to detect the target molecule. The HRP catalyzes the reaction of the DAB substrate, producing a brown color precipitate at the reaction site. This product adopts a special formula for high sensitivity, low background, good repeatability and stable storage. It is usually recommended to add 80-100uL DAB Working Solution to each tissue, and 40-50uL DAB Working Solution to each square centimeter of imprinting film.

The 60mL of Enhanced DAB Chromogenic kit (20x) can be used for staining approximately 600 tissue sections or staining imprinted membranes with an area of 1200 square centimeters.

| | |
|-------------------|---|
| Catalog No. | 56280007 |
| Size | 3mL / 10mL |
| Product Category | Immuno Detection |
| Kit Components | 1) Solution-A: 3mL/10mL 2) Solution-B: 3mL/10mL 3) Solution-C: 60mL/2 x 100mL |
| Storage/Stability | -20°C / 1 year |