



## Rapid TA Cloning kit

---

This product is based on TA cloning vector for efficient and convenient PCR cloning. The ligation reaction can be completed in 5 minutes at 25°C. The TA Cloning Vector is strictly purified, and the 3'-protruding T base is stable, and the vector self-linking probability is low. The number of blue spots for 0.5-2 kb inserts can be less than 5%, and the positive rate of white spot insertion can reach more than 90%, so blue-white spot screening can be avoided.

The kit provides DNA ligase, a reaction buffer and a 500bp PCR product (control insert) as a positive control for the ligation reaction. It is recommended to use M13 Forward and M13 Reverse or T7 Promoter and T3 Promoter universal primers for DNA sequencing.

Kit components:

1. T Cloning Vector (20μL)
2. 500bp Control Insert (5μL)
3. DNA Ligase (12μL)
4. 2x Ligation Buffer (300μL)
5. User Manual (1)

Catalog No.	515001
Size	20 preps
Product Category	Cloning / Ligation
Storage/Stability	-20°C/1 year
Shipping	Gel Packs