



Quick Western Transfer Buffer [10x] (powder)

Quick Western Transfer Buffer [1x] is an efficient, safe and non-toxic protein transfer product. The heat generated during the transfer process is less than that of ordinary transfer buffer, and no additional ice-water bath is needed. If the detection time is long, just use an ice box to cool down in the transfer tank. It can quickly complete the entire film transfer process within 5-30 minutes.

This product is safe and non-toxic, and the powder form is convenient for storage. Before use, it needs to be dissolved in pure water, and 10% absolute ethanol is added, and the final volume is adjusted to 1L. 10x solution should be prepared for storage. Before use, dilute to 1x and add absolute ethanol to a final concentration of 10%.

Catalog No.	503013
Size	1L / 10 x 1L
Product Category	Western-Blot
Storage/Stability	Ambient/1 year
Shipping	Ambient